

# Specialty PHARMACY


2104 AL Hwy 157 Cullman, AL 35058 256-734-3146

[www.specialtypharmacyonline.com](http://www.specialtypharmacyonline.com)

Read all of this information and instructions carefully before dispensing cream from the Topi-Click™.

1. Pop off the cap by placing thumb in groove on the front and pressing upward.
2. Turn the base in direction shown below in picture #2 for “2 clicks” to dispense one dose or \_\_\_ number of clicks per topical provider’s instructions.
3. Apply cream by rubbing applicator pad to area of skin as directed by topical provider.
4. Replace cap.
5. Once you’ve completed all of your prescribed doses, dispose of the Topi-Click™ – it is not refillable.

**\*The Refill Indicator** (See Fig. 5 below): If your prescription has refills, the Topi-Click™ will indicate when it’s almost empty and time for you to call the pharmacy to get your refill ready for pick-up. The refill indicator is on the Topi-Click™ cylinder near the top. When the blue or pink colored indicator inside the cylinder rises to the word **REFILL**, it’s time to call the pharmacy above for a refill.



The Topi-Click™ is a disposable dispenser and is not refillable.